

Product Name	Optilube™, sachets of lubricant		
Reference numbers	5469000		
Product Type	Optilube™ Lubricating jelly		
Manufacturer	Intersurgical		
Materials	Bleached Kraft Paper Surlyn (Ionomer)	PE Water soluble sterile lubricating jelly	Aluminium Foil
Packaging	Individually packaged		
Unit Packaging	Individually packaged		
Box quantity	150		
Pallet quantity	38400		
Instructions for Use	Instructions for use with product	NO	
	Instructions for use with each box	YES	
Certification	CE0120	Class 2A	
Applicable Standards	BS EN ISO 13485: Medical devices — Quality management systems — Requirements for regulatory purposes		
	ISO 9001: Quality Management Systems - Requirements		
	ISO 15223-1: Medical Devices - Symbols to be Used with Medical Device Labels, Labelling, & Information to be Supplied		
Sterile	Yes		
Shelf life	5 Years		
Disposal Methods	As hospital guidelines		
Re-use	 Single use only		

Product characteristics

Technical information

Weight (g)	5
Seal width at top (mm)	5 ± 1
Seal width at sides (mm)	7 ± 2
Seal width at bottom (mm)	7 ± 2



**Material
Abbreviations**

ABS	Acrylonitrile Butadiene Styrene
AL	Aluminium
HDPE	High Density Polyethylene
HIPS	High Impact Polystyrene
LDPE	Low Density Polyethylene
MBS	Methacrylate Butadiene Styrene
NRL	Natural Rubber Latex
PC	Polycarbonate
PE	Polyethylene
PP	Polypropylene
PVC	Poly Vinyl Chloride
SB	Styrene Butadiene
Si	Silicone
SS	Stainless Steel
TPE	Thermoplastic Elastomer
SBC	Styrene Butadiene Copolymer

